

RASA™

Rapid Assisted Screening Audiometry

MAICO Simplifying Pure Tone Hearing Screening Methods with RASA!

Traditional pure tone audiometers lack safeguards to prevent common screening errors related to frequency, hearing level, and ear selection. To address this, MAICO has integrated RASA into its product – streamlining the process and enabling faster, more accurate, and more accessible hearing screenings for all trained personnel. Traditional pure tone screening has just evolved with MAICO RASA.

What Makes the RASA Test Unique?

RASA integrates with digital platforms to display only essential interface elements during screening, reducing visual clutter and enhancing test accuracy.



Guided Workflow

Just tap, present, and record. The workflow guides you – tones are presented, responses marked, and the test advances automatically through frequencies and ear sides, streamlining the process for consistent, efficient results.

Protocol Simplicity

Customize protocols to fit your specific workflows, removing unnecessary steps like unused frequencies or manual level adjustments that can introduce variability. Testing stays fast, consistent, and precise.

Pass/Refer Results

Immediate Pass/Refer results – no audiogram interpretation needed. Get straightforward answers, right at the point of care.

RASA™: Built for Efficiency. Designed for Impact.

Confidence in Every Test

Built-in noise monitoring measures ambient noise levels, providing confidence in results, even in challenging environments.

Speed without Compromise

Screen faster than traditional methods – reach more people, stay on schedule, and deliver results you can trust. Take advantage of manual or automatic operating modes.

Made for Every Screener

RASA supports screening by any trained staff member – making hearing checks more accessible than ever.

Uniform Testing, Elevated Quality

Standardized protocols, clear screening results and reporting ensure every screening meets the same high standard.

RASA isn't just a new name – it's a smarter way to screen. Faster results, effortless use, and reliable outcomes, every time. That's the RASA difference.

Experience RASA with
easyTone



Global Contact:

MAICO Diagnostics GmbH

Sickingenstr. 70 - 71 · 10553 Berlin · Germany

Tel.: +49 (30) 70 71 46-50

sales@maico.biz · www.maico-diagnostics.com

US Contact:

MAICO Diagnostics US

10393 West 70th Street · Eden Prairie · MN 55344 · USA

Tel.: +1 (888) 941-4201 · Fax: +1 (952) 903-4100

info@maico-diagnostics.com · www.maico-diagnostics.com



12/2025